

## LAPAROSCOPE HD+2

Ø 5.5mm; 30°; 300 mm

autoclavable

**SKU:** 601-005-30HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)



30°

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)