

## LAPAROSCOPE HD+2

Ø 5.5mm; 30°; 300 mm

autoclavable

**SKU:** 601-005-30HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Telescopes</a>  |
| <b>SCOPE</b> | <a href="#">Laparoscopy</a> |



### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Telescopes</a>  |
| <b>SCOPE</b> | <a href="#">Laparoscopy</a> |