

LAPAROSCOPE HD+2

Ø 5.5mm; 45°; 300 mm

autoclavable

SKU: 601-005-40HD2

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



45°

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)