

## LAPAROSCOPE HD+2

Ø 10.0mm; 0°; 330 mm

autoclavable

**SKU:** 601-010-00HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)