

LAPAROSCOPE HD+4

Ø 5.5mm; 30°; 300 mm

autoclavable



SKU: 601-005-30HD4

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)