

LAPAROSCOPE HD+4

Ø 10.0mm; 30°; 330 mm

autoclavable

SKU: 601-010-30HD4

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



30°

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)