

LAPAROSCOPE HD+4

Ø 5.5mm; 0°; 300 mm

autoclavable

SKU: 601-005-00HD4

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)