

LASER GUIDE TUBE Ø 1.0 MM WITH RETRACTOR

for laser working element 630-048-00



SKU: 630-048-60

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Instruments for Laser Therapy](#)
TYPE [Laser Guide Tubes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Instruments for Laser Therapy](#)
TYPE [Laser Guide Tubes](#)