

## LASER GUIDE TUBE Ø 1.4 MM WITH RETRACTOR

for laser working element 630-048-00



**SKU:** 630-048-64

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Laser Guide Tubes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Laser Guide Tubes</a>