

LASER-WORKING ELEMENT

passive cutting; closed handle

for fibres up to Ø 1.2mm



SKU: 630-047-12

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Instruments for Laser Therapy](#)
TYPE [Working Elements](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Instruments for Laser Therapy](#)
TYPE [Working Elements](#)