

LASER WORKING ELEMENT FOR EXCHANGEABLE

for laser guide tubes; to use with

Ø 4.0mm/30° telescopes/300mm



SKU: 630-048-00

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements