

## LASER-WORKING ELEMENT, PASSIVE

for flexible laser-probes up to 1.2 mm

FL 302 mm; for 2.9 mm telescopes/ 22 Fr.

Resectoscopy-systems (630-304-01+630-309-01) Closed handle; LuerLock connector



**SKU:** 630-046-12

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>