

## LASER-WORKING ELEMENT

passive cutting; closed handle

for fibres up to Ø 0.8mm

Compatible with Resectoscopy-Systems 24/26 Fr. (630-204-06 /-26 & 630-209-06/-26) 27/28.5 Fr. (630-204-08/-28 & 630-209-08/-28) Urethrotome sheath 21 Fr. (630-400-21)



**SKU:** 630-047-08

**Categories:** [Endoscopy](#), [Urology](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Laser Therapy](#)

**TYPE** [Working Elements](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Laser Therapy](#)

**TYPE** [Working Elements](#)