

## LAWRENCE HAEMOSTATIC FORCEPS

curved; serrated

28.0cm (11.2")


**SKU:** 040-457-28

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ070
<b>KLS Martin</b>	13-583-28
<b>AUTHORS</b>	<a href="#">LAWRENCE</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ070
<b>KLS Martin</b>	13-583-28
<b>AUTHORS</b>	<a href="#">LAWRENCE</a>

## SPECIFICATIONS

[curved](#), [serrated](#)