

LEKSELL-STILLE BONE RONGEUR

jaw 5mm; double convertible

24.0cm (9.6")




SKU: 196-349-04

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Rongeurs
KLS Martin	23-907-24
AUTHORS	LEKSELL , STILLE
Medicon	31.51.05
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Rongeurs
KLS Martin	23-907-24
AUTHORS	LEKSELL , STILLE
Medicon	31.51.05

SPECIFICATIONS

[straight](#)