


## LERICHE HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

15.0cm (6")



**SKU:** 040-382-15

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH170
<b>KLS Martin</b>	13-328-15
<b>AUTHORS</b>	<a href="#">LERICHE</a>
<b>Medicon</b>	15.58.15
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH170
<b>KLS Martin</b>	13-328-15

## AUTHORS

[LERICHE](#)

## Medicon

15.58.15

## SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)