


LERICHE HAEMOSTATIC FORCEPS

straight; serrated

15.0cm (6")




SKU: 040-380-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH160
KLS Martin	13-326-15
AUTHORS	LERICHE
Medicon	15.56.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH160
KLS Martin	13-326-15

AUTHORS

[LERICHE](#)

Medicon

15.56.15

SPECIFICATIONS

[serrated](#), [straight](#)