

LICHTENBERG NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

SKU: 050-366-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM146
KLS Martin	20-676-20
AUTHORS	LICHTENBERG
Medicon	11.26.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM146
KLS Martin	20-676-20

AUTHORS

[LICHTENBERG](#)

Medicon

11.26.20

SPECIFICATIONS

[serrated](#), [straight](#)