

## LICHTENBERG NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

**SKU:** 050-366-20

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>Aesculap</b>       | BM146   |
| <b>KLS Martin</b>     | 20-676-20   |
| <b>AUTHORS</b>        | <a href="#">LICHTENBERG</a>                         |
| <b>Medicon</b>        | 11.26.20  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                |
|-------------------|--------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>           |
| <b>TYPE</b>       | <a href="#">Needle Holders</a> |
| <b>Aesculap</b>   | BM146                          |
| <b>KLS Martin</b> | 20-676-20                      |

## AUTHORS

[LICHTENBERG](#)

## Medicon

11.26.20

## SPECIFICATIONS

[serrated](#), [straight](#)