

LICHTWITZ ANTRUM TROCAR

with Luer-Lock connector

16 Gauge; 110mm

SKU: 084-172-00

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Trocars
KLS Martin	18-520-11
AUTHORS	LICHTWITZ
Medicon	09.52.10

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Trocars
KLS Martin	18-520-11
AUTHORS	LICHTWITZ
Medicon	09.52.10