

## LIGATURE ADJUSTER

straight

23.0cm (9.2")

**SKU:** 056-530-23

**Categories:** [Open Surgery](#), [Suture](#)

---

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Knot Tiers, Ligature Catcher, Loop Twisters</a>
<b>Aesculap</b>	BM950
<b>KLS Martin</b>	20-886-23
<b>Medicon</b>	13.67.23

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Knot Tiers, Ligature Catcher, Loop Twisters</a>
<b>Aesculap</b>	BM950
<b>KLS Martin</b>	20-886-23
<b>Medicon</b>	13.67.23