

LIGATURE ADJUSTER

straight

23.0cm (9.2")

SKU: 056-530-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Knot Tiers, Ligature Catcher, Loop Twisters
Aesculap	BM950
KLS Martin	20-886-23
Medicon	13.67.23

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Knot Tiers, Ligature Catcher, Loop Twisters
Aesculap	BM950
KLS Martin	20-886-23
Medicon	13.67.23