

LIGATURE NEEDLE

flexible

20.0cm (8")

SKU: 056-450-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Needles |
| Aesculap | BM861 |
| KLS Martin | 20-825-20 |
| Medicon | 13.11.50 |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Needles |
| Aesculap | BM861 |
| KLS Martin | 20-825-20 |
| Medicon | 13.11.50 |