

LIGHT CABLE CONNECTOR

Wolf/Dyonics; Clip

telescope side

Telescope-side
(light output)



SKU: 540-880-05

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors