

LIGHT CABLE CONNECTOR

Wolf/Dyonics; Snap

telescope side

Telescope-side
(light output)



SKU: 540-880-13

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | Connectors, Adaptors |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | Connectors, Adaptors |