

LIGHT CABLE CONNECTOR

Wolf

unit side

Projector-side
(light input)



SKU: 540-890-04

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors