

LIGHT GUIDE CABLE CONNECTOR

Olympus; Switch

light source side (new model)

Projector-side
(light input)



SKU: 540-890-14

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Connectors, Adaptors