

## LIGHT GUIDE CABLE CONNECTOR

Olympus/ACMI - Standard

light source side (old model)

Projector-side  
(light input)



**SKU:** 540-890-10

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Connectors, Adaptors</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Connectors, Adaptors</a>