

LIGHTLINE TROCAR SHORT, PYRAMIDAL,

permeable, Ø 5.5 mm, for WL 75 - 90 mm

with hollow grinding



SKU: BLUO/04-00012-S

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)