

LIGHTLINE TROCAR SLEEVE PLASTIC Ø 11 MM

maintenance free stopcock, WL 130mm

threaded tube, distal bevelled



SKU: BLUO/04-00412-1-F

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)