

LIGHTLINE TROCAR SLEEVE PLASTIC Ø 5.5 MM

maintenance free stopcock, WL90mm

threaded tube, distal bevelled



SKU: BLUO/04-00414-1-SF

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)