

## LIGHTLINE TROCAR SLEEVE PLASTIC Ø 5.5 MM

maintenance free stopcock, WL 120mm

threaded tube, distal bevelled



**SKU:** BLUO/04-00414-1-F

**Categories:** [Blue endo](#), [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)