

LIGHTLINETROCAR -SLEEVE, STAINLESS STEEL

Ø 11mm, maintenance free stopcock

WL120mm, rough tube, distal bevelling



SKU: BLUO/04-00413-2-F

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)