

## LIGHTLINETROCAR -SLEEVE, STAINLESS STEEL

Ø 11mm, maintenance free stopcock

WL120mm, rough tube, distal bevelling



**SKU:** BLUO/04-00413-2-F

**Categories:** [Blue endo](#), [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)