

LOOP ELECTRODE

angled; 24 / 26Fr.

blue; for 12° telescopes; monopolar



SKU: 630-240-05

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes