

LOUTE WIRE TIGHTENER

21.0cm (8.4")

SKU: 190-780-21

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Twisting Instruments
Aesculap	KY220
KLS Martin	22-542-21
AUTHORS	LOUTE
Medicon	76.05.50

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Twisting Instruments
Aesculap	KY220
KLS Martin	22-542-21
AUTHORS	LOUTE
Medicon	76.05.50