

M.G.HOSPITAL NEEDLE HOLDER

straight; serrated; grooved jaw

17.0cm (6.8")



SKU: 050-072-17

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM257
KLS Martin	20-348-17
AUTHORS	M.G.HOSPITAL
Medicon	10.14.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM257
KLS Martin	20-348-17

AUTHORS

[M.G.HOSPITAL](#)

Medicon

10.14.18

SPECIFICATIONS

[serrated](#), [straight](#)