

MASSON NEEDLE HOLDER

straight; serrated; grooved jaw

27.0cm (10.8")

SKU: 050-150-27

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM275
KLS Martin	20-408-27
AUTHORS	MASSON
Medicon	10.40.27
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM275
KLS Martin	20-408-27

AUTHORS

[MASSON](#)

Medicon

10.40.27

SPECIFICATIONS

[serrated](#), [straight](#)