


## MASSON NEEDLE HOLDER

straight; serrated; grooved jaw

27.0cm (10.8")

**SKU:** 050-150-27

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM275
<b>KLS Martin</b>	20-408-27
<b>AUTHORS</b>	<a href="#">MASSON</a>
<b>Medicon</b>	10.40.27
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM275
<b>KLS Martin</b>	20-408-27

**AUTHORS**

[MASSON](#)

**Medicon**

10.40.27

**SPECIFICATIONS**

[serrated](#), [straight](#)