

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

16.0cm (6.4")



SKU: 050-080-16

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM235
KLS Martin	20-390-16
AUTHORS	HEGAR , MAYO
Medicon	10.18.15
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM235
KLS Martin	20-390-16

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.15

SPECIFICATIONS

[serrated](#), [straight](#)