

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

18.0cm (7.2")



SKU: 050-080-18

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM236
KLS Martin	20-390-18
AUTHORS	HEGAR , MAYO
Medicon	10.18.18
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM236
KLS Martin	20-390-18

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.18

SPECIFICATIONS

[serrated](#), [straight](#)