

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

16.0cm (6.4")



SKU: 050-090-16

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM235
KLS Martin	20-396-16
AUTHORS	HEGAR , MAYO
Medicon	10.18.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM235
KLS Martin	20-396-16

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.15

SPECIFICATIONS

[serrated](#), [straight](#)