

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

18.0cm (7.2")




SKU: 050-090-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM236
KLS Martin	20-396-18
AUTHORS	HEGAR , MAYO
Medicon	10.18.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM236
KLS Martin	20-396-18

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.18

SPECIFICATIONS

[serrated](#), [straight](#)