

MAYO-HEGAR NEEDLE HOLDER

straight; serrated; grooved jaw

14.0cm (5.6")



SKU: 050-080-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM245
KLS Martin	20-390-14
AUTHORS	HEGAR , MAYO
Medicon	10.14.14
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM245
KLS Martin	20-390-14

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.14.14

SPECIFICATIONS

[serrated](#), [straight](#)