

MC PHERSON BIPOLAR COAGULATION FORCEPS

straight; blunt; Ø 0.45mm

10.0cm (4"); flat plug



SKU: 022-210-45

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	blunt , straight