

MC PHERSON BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.25mm

10.0cm (4"); flat plug



SKU: 022-211-25

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	angled , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	angled , sharp