

MC PHERSON BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 0.45mm

10.0cm (4"); flat plug

SKU: 022-211-45

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
AUTHORS	MC PHERSON
SPECIFICATIONS	angled , blunt