

## MC PHERSON BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 0.45mm

10.0cm (4"); flat plug

**SKU:** 022-211-45

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>AUTHORS</b>	<a href="#">MC PHERSON</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>AUTHORS</b>	<a href="#">MC PHERSON</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>