

## MC PHERSON BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.25mm

10.0cm (4"); flat plug



**SKU:** 022-210-25

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>AUTHORS</b>	<a href="#">MC PHERSON</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>AUTHORS</b>	<a href="#">MC PHERSON</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>