

MENGHINI PUNCTURE NEEDLE

Ø 1.8 x 168mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-117-18

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
KLS Martin	18-710-18
Medicon	09.75.78

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
KLS Martin	18-710-18
Medicon	09.75.78